



# ReaScan® Revolution

# Welcome to the ReaScan Revolution

Proven lateral-flow technology and proprietary assay development techniques are integrated into the next generation ReaScan+ product family. Its ease of use combined with result traceability and LIS connectivity add confidence to the test results.





#### The ReaScan Revolution

ReaScan+ reader is a compact rapid test device with an intuitive graphical user interface. It combines rapid and reliable results with ease of use. ReaScan+ reader offers the option of Immediate Reading or Incubation Mode (walk away) allowing flexible workflow in the laboratory.

Up to 300 results are stored into the internal non-volatile memory

Touch screen with an easy to use User Interface supporting multiple languages

The
2D code on
ReaScan cassettes
is automatically
scanned to verify the
correctness of the
test



User and
Patient ID, date
and time and test
kit batch number are
stored with each
result for maximum
traceability

Results can be automatically transferred to the LIS or laboratory middleware

Optional 1D/2D barcode scanner to automate the input of sample or patient ID

# ReaScan+ C6 LYME IgG

Reagena is a pioneer in rapid diagnostics for zoonotic diseases and lately focused on the diagnostic tests for Lyme Borreliosis and Lyme neuroborreliosis. ReaScan+ C6 LYME IgG is a lateral flow rapid test intended for detection of Borrelia burgdorferi specific IgG antibodies in human serum and cerebrospinal fluid (CSF). The test is intended for use as an aid for diagnosis of Lyme Disease.

# ReaScan+ LYME IgM/IgG Combo

ReaScan+ LYME IgM/IgG Combo is a lateral flow rapid test intended for qualitative detection of Borrelia burg-dorferi specific antibodies in human serum and cerebrospinal fluid (CSF). The test uses C6 (VIsE) and PepC10 (OspC) antigens in separated test lines. PepC10 generates an early IgM response and thus it is likely to react especially in acute Lyme disease.

Product	Test performance	Packaging	Ref. number
ReaScan+ C6 LYME IgG	Specificity 99 % (CSF), 95 % (serum) Sensitivity 97 % (CSF), 97 % (serum)	10 tests	114313
ReaScan+ LYME IgM/IgG Combo	PPA 94,8 % (CSF, serum/plasma) NPA 95,6 % (CSF, serum/plasma)	10 tests	115103

The diagnostic sensitivity of the ReaScan+ C6 LYME IgG was determined by analyzing 108 paired serum and CSF samples from definite LNB patients. The diagnostic specificity of the ReaScan+ C6 LYME IgG was evaluated by analyzing paired serum and CSF samples from 104 non-LNB patients. EQV results counted as positive.

ReaScan+ LYME IgM/IgG Combo PPA and NPA values for 326 samples including 168 positive serum/plasma samples, 117 negative serum/plasma samples, 23 positive CSF samples and 18 negative CSF samples. EQV ReaScan+ LYME IgM/IgG Combo results were calculated as negatives.

#### ReaScan CXCL13

ReaScan CXCL13 is a lateral flow rapid test for semiquantitative determination of CXCL13 in human cerebrospinal fluid CSF, which helps clinicians in the treatment decision for suspected Lyme neuroborreliosis (LNB) patients.

Product	Test performance	Packaging	Ref. number
ReaScan CXCL13	Specificity 96 % Sensitivity 100 % with cut-off 250 pg/mL	10 tests	114253

Performance was obtained by comparing ReaScan CXCL13 with a widely used CXCL13 ELISA (Human CXCL13/BLC/BCA-1 Quantikine, R&D Systems, Minneapolis, USA) with 225 patient samples with a suspected diagnosis of LNB. Ref.: Pietikäinen A, et al, Point-of-care testing for CXCL13 in Lyme neuroborreliosis, Diagn Microbiol Infect Dis (2018), https://doi.org/10.1016/j.diagmicrobio.2018.02.013.

# ReaScan TBE IgM

ReaScan TBE IgM is a lateral flow rapid test for qualitative detection of tick-borne encephalitis (TBE) virus specific IgM antibodies in human serum and cerebrospinal fluid (CSF).

Product	Test performance	Packaging	Ref. number
ReaScan TBE IgM	Specificity 100 % (CSF), 98 % (serum) Sensitivity 92 % (CSF), 99 % (serum)	10 tests	114106

Performance was evaluated with 172 serum samples from earlier diagnosed TBE patients and 141 TBEV IgM negative serum samples. See Albinsson B et al. 2020. Multi-laboratory evaluation of ReaScan TBE IgM rapid test, 2016 to 2017. Euro Surveill. 2020 Mar;25(12).

#### ReaScan Hantavirus

Reagena is an undisputed leader in Hantavirus rapid diagnostics. Reagena offers rapid tests intended for qualitative detection for Puumala virus and Dobrava or Hantaan virus specific IgM antibodies from human serum or plasma. Test results are used as an aid in the diagnosis of acute Hantavirus infections.

Product	Test performance	Packaging	Ref. number
ReaScan+ PUUMALA IgM	Specificity 99 % Sensitivity 98 %	10 tests	115105
ReaScan DOBRAVA-HANTAAN IgM	Specificity 97 % Sensitivity 95 %	10 tests	114102

Performance of each test was obtained when evaluating test performance with panels of confirmed patient samples.











#### **Technical data**

Display: 4.3" Touch screen

Memory: 300 data sets

Connectivity: 3x USB ports, 1x LAN port

Power Supply: 12 Vdc (adapter included)

Size: 140 x 180 x 150 mm

Weight: 1150 g

Operating temperature: +15°C ... +35°C

One year limited warranty



### Feature highlights

- ♦ Fast, accurate and objective results
- **▲** Traceability of test results
- ▲ Robust and maintenance free device
- Open and flexible connection to different LIS systems
- **▲** Intuitive Graphical User Interface
- Immediate reading mode or delayed reading after incubation time

#### **Order information**

Product name	Ref. number
ReaScan°+ reader	105100
ReaScan®+ connect server	105130
Zebra DS2208 SR USB-Kit Black barcode scanner	105110
Dymo LabelWriter 450 printer	105120







